## Livento® invest

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N° de commande 083 739

## Livento® invest

Special investment material for lithium disilicate and other pressable ceramics. Fabrication of refractory stumps. Metal casting technique.

### Indications

For rapid and conventional preheating techniques.

- Investment material for lithium disilicate and other pressable ceramics.
- Fabrication of refractory stumps.
- Metal casting technique\*.

\*Note: Basically possible. Depending on the shape and size of the inner geometry of the crown, the fit of the crown on the stump may be too tight, and post-processing may be necessary.

### Preparing the wax patterns

Once the sprues have been attached, apply wetting agent and allow it to evaporate or rinse off the residue with water (approx. 30°C).

### Preparing the casting cylinder

Select the correct size of casting cylinder. Mask base and inner surface of the ring with silicone spray for easy removal of the pressing cylinder.

### Mixing ratio

Mixing ratio – 26 ml liquid: 100 g powder

### Recommended mixing liquid concentrations

These figures are to be considered guidelines only. The results can vary slightly, depending on the batch, the type of pressing ceramic and the method of work.

### Concentrate mixtures

- |         |   |
|---------|---|
| 40–60 % | Expansion fluid for inlays, onlays, depending on preparation and size |
| 60–70 % | Expansion fluid for single crowns                                     |
| 70–85 % | Expansion fluid for posterior and anterior tooth bridges              |
| 100 %   | Expansion fluid for refractory stumps                                 |
| 100 %   | Expansion fluid for metal casting technique                           |

**Please note:** The higher the amount of concentrate in the liquid, the higher the expansion values.

### Investment ring loading

- |       |                          |
|-------|--------------------------|
| 100 g | max. 3 similar wax parts |
| 200 g | max. 5 similar wax parts |
| 300 g | max. 7 similar wax parts |

### Mixing

- |        |   |
|--------|---|
| Step 1 | Pre-fill with fluid mixture   |
| Step 2 | Sprinkle with powder  |
| Step 3 | Allow to soak in for approx. 10 seconds   |
| Step 4 | Mix vigorously for 10–15 seconds  |
| Step 5 | Stir for 1 minute under vacuum at 300–350 rpm                                   |
| Step 6 | Maintain vacuum for a further 15 seconds  |
| Step 7 | Topping up the press investment ring is performed at the lowest shaker setting. |
- Once the investment ring has been filled, shaking is stopped.

Processing time approx. 5 min.

Processing temperature 19–22°C

**Please note:** Room and storage temperatures exceeding 22°C curtail the working time! Avoid creating dust. Do not inhale dust (wear a face mask).

### Burnout/Preheating

#### a) Rapid preheating:

After investment, the 100 g cylinders should be allowed to harden at room temperature for at least 23–28 minutes, the 200 g and 300 g cylinders for at least 25–30 minutes.

Note: When using models from 3D printing or milling systems, a longer hardening time is recommended. Allow the 100 g cylinder to harden at room temperature for 35–40 minutes and the 200 g and 300 g cylinders for 40–45 min. After 20 minutes carefully remove the lid of the investment ring and the silicone ring (leave on the base). After removing the acrylic base, place the cylinder in the preheated furnace at 850°C.

#### Attention:

Do not open the preheating furnace during the heating phase, as the generated wax vapours can ignite due to draught.

#### b) Conventional preheating (overnight preheating):

As the full setting expansion is achieved with the overnight technique, the fluid concentration should be reduced by 5–10 %, depending on the ceramic and the work.

Recommended heating rate: 3–5°C/minute for linear furnace control.

If several investment rings are placed in the furnace, the dwell time at the final temperature of 850°C should be extended by 10 minutes per additional investment ring.

#### Attention:

- In case of a ribbed floor (furnace without floor heating!), the investment ring can be placed directly with the opening facing downwards.
- In case of flat floors, please make sure that the wax burn-out occurs outside of the investment ring, e.g. by tipping the investment ring in direction of the rear wall.
- In case of furnaces with floor heating, please ensure that the investment ring is placed in a distance of approx. 10 mm from the floor.

Preheating and dwell times at 850°C

Max. thermal stability: 1200°C

100g cylinder	at least 45 minutes
200g cylinder	at least 60 minutes
300g cylinder	at least 70 minutes

Then the press cylinder is filled with the desired pressing ceramic in the usual way and placed in the pressing furnace.

#### Cooling and demolding of the cylinder

After cooling to room temperature, the cylinder will be gently removed from the mold and then sandblasted with glass beads and a pressure of 2.5–3.0 bar.

#### Fabrication of refractory stumps

1. Mixing the investment material with 100% expansion liquid.
2. Demould the stumps after 30 min. setting time of the investment material.
3. Perform degassing firing in the ceramic furnace.

#### Program:

- Drying: 5 min.
- Closing time: 5 min.
- Preheating temperature: 180°C
- Temperature 1: 700°C
- Temperature increase: 10°C per min.
- Hold time: 10 min. (without vacuum)
- Temperature 2: 1,050°C
- Temperature increase: 55°C per min.
- Hold time: 5 min. (without vacuum)

#### Physical properties (pure mixing liquid)

Total linear expansion: 3%

#### Important

These instructions supersede all previous editions.

#### Danger

H373 May cause damage to organs through prolonged or repeated exposure.

#### Warning

P260 Do not breathe dust/fume/gas/mist/vapours/spray.  
 P284 Wear respiratory protection.  
 P314 Get Medical advice/attention if you feel unwell.

**Product contains:** quartz and cristobalite.

For further information ask for the Health and Safety Datasheet from Cendres+Métaux SA.

**Storage:** Livento® invest cannot be stored for an unlimited period, please refer to the back of the bag. As investment compounds are sensitive products, they must be protected against air, moisture and warmth in their original packaging. The ideal storage temperature is 18°C.

 Attention: store at room temperature!

The **special liquid** is sensitive to low temperatures such as when being transported during cold weather.

**Traceability, lot numbers:** When placing complaints, always quote the lot number (back of bag).

**Allergies:** Certain constituents of investment compounds may cause contact allergies on skin while using the investment. It is therefore advisable to wear gloves to protect the hands, and a face mask.  
 Never mix investment compounds with different lot numbers.

#### Liability

The information provided herein represents state-of-the-art technology and assumes that the product has been checked carefully. We guarantee impeccable quality but are not liable for results due to further processing which is not usually within our control.

The manufacturer assumes no responsibility whatsoever for damage due to non-adherence to these instructions, especially if the investment compound is used after the expiry date.

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