

Instructions for use

Description

Grinding points based on ceramic bonded refined corundum.

Indication

Efficient high quality cutting abrasives, used for casting alloys for porcelain bonding alloys, containing more than 50 % gold.

Characteristics

- refined corundum, ceramic bonded
- excellent cutting efficiency
- fine cutting grain structure
- wide range of cutting speed (5000–12000 rpm)
- true running due to precision centred shaft (∅ 2.35mm)
- long use and constant grain structure during wear

Instructions for use

Work out the metal surface only in one direction with moderate pressure.

Safety precautions

Do not breathe dust

Presentation
box with 12 points
08052192 Cerasiv brown

Rx only