

Portfolio.

Order No.	Product	Description	Kit
08056952	Livento® press – Soprano® 10	Bleaching Extended Kit	
08056950	Livento® press	Bleach ET	3 pc.
08056942	Livento® press	Bleach 1 MT	3 pc.
08056943	Livento® press	Bleach 2 MT	3 pc.
08056944	Livento® press	Bleach 3 MT	3 pc.
08056945	Livento® press	Bleach 4 MT	3 pc.
08056946	Livento® press	Bleach 1 LT	3 pc.
08056947	Livento® press	Bleach 2 LT	3 pc.
08056948	Livento® press	Bleach 3 LT	3 pc.
08056949	Livento® press	Bleach 4 LT	3 pc.
08056954	Livento® press	Bleach MO 0	3 pc.
08056955	Livento® press	Bleach MO 1	3 pc.
08056951	Livento® press	Bleach HO 0	3 pc.
08056953	Soprano® 10	Bleach Opaque Dentine	15 g
08056956	Soprano® 10	Bleach Dentine 1	25 g
08056957	Soprano® 10	Bleach Dentine 2	25 g
08056958	Soprano® 10	Bleach Dentine 3	25 g
08056959	Soprano® 10	Bleach Dentine 4	25 g
08056960	Soprano® 10	Cuspid	15 g
08056961	Soprano® 10	Bleach Enamel 1	25 g
08056962	Soprano® 10	Bleach Enamel 2	25 g
08056963	Soprano® 10	Enamel EO	25 g
08055214	Soprano® 10	Enamel E1	25 g
08055218	Soprano® 10	Enamel Clear	25 g
08056964	Soprano® 10	Bleach Transpa-white	15 g
08055179	Soprano® 10	Enamel Effect light-blue	15 g
08055219	Soprano® 10	Opal Intensive 2	15 g
08055178	Soprano® 10	Mamelon cream	15 g
08055277	Soprano®	Modeling Liquid	100 ml
08055283	Soprano®	Glaze Fluid	25 ml
08055184	Soprano®	Glaze Paste	5 g

esthetic.line

Livento® press

Soprano® 10

Bleaching Extended Kit.

New
&
unique



The new Bleaching Extended Kit

Part of the lithium disilicate press ceramic system Livovento® press and Soprano® 10

Your palette of Livovento® press and Soprano® 10 has just got better – As a result of high market demand from Middle East and American markets using more than 70% brighter than B1, the Bleaching Extended Kit is here!

Our key opinion leader teams globally and especially in the Middle East just simply needed a great range of colours. To enhance restorations from monolithic crowns to copings with multi layering, high end aesthetics, Soprano® 10 is now reaching even higher levels in aesthetically demanding cases.



The Bleaching Extended Kit contains:

- **Livovento® press Bleach LT**
low transparency in four bleach shades, from high value Bleach 1 to high chroma Bleach 4.
- **Livovento® press Bleach MT**
medium transparency in four bleach shades, from high value Bleach 1 to high chroma Bleach 4.
- **Livovento® press Bleach MO**
medium opacity ingot available in two shades 0 and 1, used in some special cases that need to mask more than the LT.
- **Livovento® press Bleach HO**
high value, high opacity, usually used to mask the underneath stub like metal abutment or very dark tooth colour.
- **Livovento® press Bleach ET**
Enamel Bleach like ingots, useful for prep-less and minimal prepared cases, keeping the original teeth colour or slightly brighter.
- **Soprano® 10 Bleach Dentine**
we covered all bleaches dentine colours from 1 to 4, that will help in layering technique with ease.
- **Soprano® 10 Bleach Enamel**
special enamel for bleaches.
- **Soprano® 10 Enamel and Effects**
a group of enamel and enamel effects works perfectly with bleaches giving a natural look.
- **Soprano® Glaze**
a natural glaze paste with liquid.